



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
DIAMOND WHEELS

**Product Code :**  
NMCBS139040



**Description :**

DIAMOND WHEELS For long life reverse / face direction of the wheel after 1 to 2 month The hole / bore of the wheel should fit in bearing position of the shaft Regular use of proper dressing stick on wheel increases cutting rate and also increases life of wheel DIAMOND WHEELS, ENT Equipment exporters, ENT Equipment manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)