

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Male External Catheter

Product Code :
NMCBS35034



Description :

Male External Catheter Position Over Tip of Penis Then Unroll the latex sleeve over entire length. Use adhesive strap around top of sleeve to secure to penis avoiding constriction. Afix distal end to drainage tubing. Used for Urine incontinence in male patients. Made from natural latex rubber. Sterile, Individually packed in peelable pouch pack

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)