

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fluoride Oxalate Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
NMCBS35090



Description :

Fluoride Oxalate Tubes (MOQ: 100,000 Pcs Assorted Size) This PET Fluoride/Oxalate Glucose Blood Collection tube is utilized for glucose determinations. Glucose values in unpreserved blood samples deteriorate swiftly after collection as glucose is used by the blood cells. The Fluoride will stop enzymatic activity while oxalate acts as an anti-coagulant. Mixing recommendation: Fluoride/Oxalate tubes should be gently inverted 180° and back 8-10 times. Centrifugation instructions: 1300g for 10 minutes at 25° C room temperature

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)