

Product Name :
Autoclave Steam Indicator Tape**Product Code :**
NMALV10037**Description :**

Autoclave Steam Indicator Tape Light Yellow indicator tapes has stripes that changes colour from yellow to dark brown/black when passes through steam sterilization Cycle (i.e. 20 minutes of exposure to a temperature of 121½ c or 5 minutes of exposure to a temperature of 134½C)

Technical Specification :

Autoclave tape is used as a process indicator in steam sterilizers. It is very well suited for the clinical and research environments. It features an aggressive, pressure-sensitive adhesive designed to adhere to a variety of wraps so that the pack stays secure during sterilization and can be used on direct contact with substrate. Tape also features a chemical indicator to verify that it has gone through the sterilization process.