

Email: sales@naugramedical.com

Phone: 9896600003

Product Name:

Thermal Conductivity Of Metal Rod

Product Code: NMCBS11160



Description:

Thermal Conductivity Of Metal Rod The experimental set up consists of metal bar, one end of which is heated by an electric heater while the other end of the bar projects inside the cooling water jacket. A cylindrical shell filled with the asbestos insulating powder surrounds the middle portion of the bar. The temperature of the bar is measured at different sections. Heat Input to the heater is given through variac and measured by Digital Voltmeter & Digital Ammeter. By varying the heat input rates, wide range of experiments can be performed. Water under constant head condition is circulated through the jacket and its flow rate and temperature rise is noted. Experiments To plot the temperature distribution along the length of Bar To determine the thermal conductivity of given bar at various temperatures Utilities Required Water supply 3 lit/min (approx.) Drain Electricity Supply: 1 Phase, 220 V AC 2 Amp. Table for set-up support Thermal Conductivity Of Metal Rod, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

1/1