

**Product Name :**  
Binocular Stereo Microscope AEI-A**Product Code :**  
NMCBS130002**Description :**

Binocular Stereo Microscope AEI-A The Binocular Stereo Microscope is designed for younger students of elementary education. It meets the needs of Science teaching and demonstration to beginners at school level. Being light weight & handy, it helps carrying around and is simple to use, yet having features of a compound Microscope Precision Optical Glass Lenses and metal components are engineered for years of reliable performance. Observation Head - Binocular head inclined at 45-deg. with adjustable inter-pupillary distance. Focussing - By Rack and Pinion mechanism. Base - Fitted on transillumination base suitable for artificial and transmitted arrangement and elbow hand rest. Magnification - Standard magnification 20X & 40X or 10x & 30x. Objectives - 2 sets of objectives 2X and 4X, selected by sliding nosepiece. Eyepieces - Standard extra wide field eyepieces WF-10X or 15X are provided. Binocular Stereo Microscope AEI-A, Binocular Loupes exporters, Binocular Loupes manufacturers

**Technical Specification :**

Industrial Stereo Microscope is a kind of binocular stereo microscopes, which can magnify micro objects continually and show clear stereo up-right picture. It can be used for observation studies in medical and health, farming and forestry, as well as public security departments, school, and scientific research institutes, and is also used for inspections, assembling and repair of tiny precision spare parts in electronics and precision machine industries.

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**