

**Product Name :**  
Video Colposcope**Product Code :**  
NMCBS134002**Description :**

Video Colposcope Image Processor : 1/6" Micro lens Color Digital CCD Pixels : 10,000 Resolution : >750 CCD Lines Electronic Shutter : Automatic S/N Ratio : >48 DB Illumination : High MCD Super Bright White Shadows Light Color Temperature : >7000° K Test Timer : Yes Up & Down Movement : Yes Focus : Automatic or Manual Focus Distance : 15 ~ 500mm Zoom : 1 ~ 40X Viewing Area : 170mm ~10mm Depth Of View : 200mm ~ 5mm Magnification Display : Yes Video Output : BNC, SVHS and DV Output (Optional) Safety Standard : Complies to IEC 60601 - 1 Gamma Processor : YES - For Operative Functions Remote Access : For Technical Services Power Supply : 100-240 VAC, 50/60 HZ Video Colposcope, Gynecology Equipment exporters, Gynecology Equipment manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,