

Product Name :
Fiber Optic Otoscope

Product Code :
NMCBS134004



Description :

Fiber Optic Otoscope Otoscope with Fiber Optic (F.O.) illumination constructed of light, high-impact Polycarbonate and 3x magnification scratch resistant acrylic viewing window. Output for pneumatic membrane test. Otoscope Set complete with 10 mini Speculum Size 2.5mm and 4mm. Set complete with C size Dry battery Handle. Fiber Optic Otoscope, Gynecology Equipment exporters, Gynecology Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)