



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Keratometer

Product Code :
NMCBS137010



Description :

Keratometer Measuring System : External reading Ocular : 15x-high-eye-point optical system Measuring range: Radius of curvature 9.4 ~ 6.4mm (in 0.05 mm steps) Diopter 36-52D(in 0.25D steps) Astigmatic axis: 0o180o(in 5o steps) Focusing distance : 76mm Coincidence focusing system : Positive measurements to .02mm accuracy Dimensions weight : 16"(h) * 7"(w) * 19"(D) and weight 7 KGS Keratometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,