

**Product Name :**  
Keratometer

**Product Code :**  
NMCBS137010



**Description :**

Keratometer Measuring System : External reading Ocular : 15x-high-eye-point optical system Measuring range: Radius of curvature 9.4 ~ 6.4mm (in 0.05 mm steps) Diopter 36-52D(in 0.25D steps) Astigmatic axis: oo180o(in 50 steps) Focusing distance : 76mm Coincidence focusing system : Positive measurements to .02mm accuracy Dimensions weight : 16"(h) \* 7"(w) \* 19"(D) and weight 7 KGS Keratometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers