

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Halogen Headlight

Product Code :
NMCBS137016



Description :

Halogen Headlight Halogen Head light for true tissue color and consistent. Long lasting Illumination. Iris diaphragm regulates diameter of Light Spot. The Halogen illumination of the lamp offers an optimal light intensity which may be regulated by a switch. The practical plastic design offers easy handling due to low weight. Head may be exchanged without removing the headband. The plug transformer 6v supplies the required power. Halogen Headlight, Ophthalmic Equipment, Ophthalmic Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)