

Product Name :
Indirect Ophthalmoscope**Product Code :**
NMCBS137017**Description :**

Indirect Ophthalmoscope Halogen Binocular Indirect Ophthalmoscope is an excellent instrument for all standard applications, as well as with uncooperative pediatric patients & in cases where the papillary diameter is less than 4mm. The need for pupillary dilatation is virtually obviated when using the SPF-capability of this dual purpose clinical instrument. Indirect Ophthalmoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,