

**Product Name :**  
Infusion Pump**Product Code :**  
NMCBS137018**Description :**

Infusion Pump - Accept injectors of 10ml, 20ml, and 50ml, identify and display the injector's type automatically, press keys to choose different brands of injectors. - Accurate and sensitive pressure sensor which shortens alarm's response time and prevents fast injection due to excessive pressure. - Unique injection supervision function. - Excellent injection operation, even small - User-friendly, easy to operate, safe, convenient, English indicators, easy to identify Infusion Pump, Ophthalmic Equipment, Ophthalmic Equipment manufacturers