

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Otoscope

Product Code :
NMCBS137025



Description :

Otoscope 2.5v Halogen Otoscope with Metal Head in Mat finishes with Pin contact fitting. Optical glass viewing window with 3x magnification. Movable Swiveling viewing window, permanently-attached. Metal reusable speculum size 2, 3, 4mm. Output port for pneumatic test. Rheostat to adjust the intensity of light. C size Battery Handle. Available in Wooden Case or in zipper pouch Otoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)