



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Slit Lamp

**Product Code :**  
NMCBS137027



**Description :**

Slit Lamp Our slit lamp is remarkable for exceptional optical performance. We make them available in four outstanding models that are quite appropriate in different kinds of medical and scientific experiments. Our slit lamp biomicroscope is available with high quality optical glass. These slit lamp microscopes are extraordinary for their clarity and detailed view of the requisite samples to be examined. They are accessible with XYZ joystick that gives convenient control for their effective performance. And when it comes to the leading slit lamp manufacturers, you can always count on our services. Further, we are also one of the key slit lamp exporters from India. Slit Lamps utilize high quality Optic, providing clarity and detail that will surprise even the most discriminating operator. The XYZ joysticks, convenient controls, and ergonomic design, combined with top quality materials and meticulous assembly, ensure year of comfortable, enjoyable use. The outstanding optical Performance, features and quality make the Kashmir Surgical series comparable to the world's elite lamps Slit Lamp, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,