
NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Riester Schiotz C Tonometer

Product Code :
NMCBS137041



Description :

Riester Schiotz C Tonometer High quality agate bearing for an extremely long service life Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions Perfect reading of the scale with red pointer All vital parts are made of stainless steel, other parts are chrome-plated Complete in a black deluxe case with velvet-look inserts Riester Schiotz C Tonometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)