

Product Name :
Automatic Capsule Loading Machine**Product Code :**
NMCBS138004**Description :**

Automatic Capsule Loading Machine We have with a range of automatic capsule loading machines that are extensively used in combination with manually operated capsule filling machines. In this machine, the capsules are automatically loaded in loading tray in cap up - body down position. Our range designed to meet the application requirements of medium and large scale pharmaceutical companies. Precision engineered, these machines are easy to operate, maintain, durable and are highly cost effective. Further, we stringently checked these on parameters in order to meet the international quality standards. Features : Easy to operate Sturdy in construction High performance Low maintenance Low noise Automatic Capsule Loading Machine, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#),