

**Product Name :**  
VIBRO SIFTER. MODEL MVS-30 (GMP MODEL)**Product Code :**  
NMCBS138010**Description :**

VIBRO SIFTER. MODEL MVS-30 (GMP MODEL) The unit is covered with S.S. 304/316 q. The unit is supplied complete with any one screen up to 60 mesh size, without dust cover, self balancing Vibro Screen producing adjustable multiplane gyratory motion, designed for floor mounting. The Vibrating Screen assembly is isolated by a set of special designed springs and is complete with specially designed vibratory motor suitable for operation on 415/440 V. , 3 Phase, 50 Cycles, AC Power supply with adjustable top, and bottom eccentric weights to get required vibrations, having all contact parts in S.S. 304 q. driven by 0.5 HP motor, size 30" dia. All M.G. visible parts are made out of S.S. or clad with S.S. The unit will be as per GMP Standard.. VIBRO SIFTER. MODEL MVS-30 (GMP MODEL), Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,