

Product Name :
TABLET DISINTEGRATION TEST MACHINE USE
STANDARD**Product Code :**
NMCBS138017**Description :**

TABLET DISINTEGRATION TEST MACHINE USE STANDARD (As per USP Standard) suitable to take two tests simultaneously with two hot plates and 2 sets of basket assembly complete with 12 guided disc TABLET DISINTEGRATION TEST MACHINE USE STANDARD, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

Technical Specification :

An instrument is equipped with the single stepper motor.
Five programmable Tests are provided to Store DT time and Temperature.
Digital indicator for temperature with optically isolated electronic heater driver with the water bath.
The Instrument is also equipped with double channels Temperature system to ensure precise temperature control of the Tank and Beakers.
Illuminator for clear observations Of Disintegration process has been provided.
Soft membrane keypad and splash proof panel.
A four line informative LCD display is provided to observe temperature and set time at a glance.
The facility provided to prevent over the shooting of temperature.

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]

Medical Hospital Equipments Manufacturers