

**Product Name :**  
RECATION VESSEL**Product Code :**  
NMCBS138022**Description :**

RECATION VESSEL Stainless Steel 304/316L quality.  $\frac{1}{2}$  MS/SS Jacket or limpet coil for steam heating / cooling.  $\frac{1}{2}$  Reaction Vessel can be provided with different types of Agitator like Anchor type / Propeller type / Paddle type Pitched blade turbine type of with high speed homogenizer.  $\frac{1}{2}$  Shaft entry sealed by Oil seal / Plain or Water cooled Stuffing Box / Single or Double Mechanical seal.  $\frac{1}{2}$  High density glass wool insulation & S.S. Cladding can be done on demand.. RECATION VESSEL, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

**Technical Specification :**

Type: Reaction Vessels  
Material: Stainless Steel  
Shape: Cylinder Shape, Round  
Color: Brown, Golden-Grey, Grey, Light Green, Metallic, Shiny Silver, Silver  
Surface Finishing: Chemical Coated, Powder Coated  
Feature: Anti Corrosive, Durable, Eco-Friendly, High Quality, Shiny Look, Stable Performance

