

**Product Name :**  
BOTTLE WASHING AND BRUSHING MACHINE**Product Code :**  
NMCBS138025**Description :**

BOTTLE WASHING AND BRUSHING MACHINE The unit is supplied complete with stainless steel soaking tank size 75 x 50 x 30 cms, 12 Nos. stainless steel rinsing jets and twin head bottle brushing unit driven by 1½ HP Motor (Single Phase A.C. 230 volts). BOTTLE WASHING AND BRUSHING MACHINE, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

**Technical Specification :**

Bottle washing machine is made of sturdy mild steel sheet and covered with thick layer of epoxy based powder coating paint.

it is having variable control upto 1500rpm.

A drill chuck is fitted for the alignment.

A special acrylic cover and enamel tray is provide for protection. It can wash upto 2000 ml bottle.

