

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Endo Camera (Digital Video Camera)

Product Code :
NMCBS139006



Description :

Endo Camera (Digital Video Camera) It is very small and have a high source of light. It is highly flexible and has zooming features. It is lightweight and economic. Specifications Pick Up Device : 1/4" CCD Pixels : PAL = 522(H) X 582(V) Horizontal Resolution : NTSC/PAL = 450 Lines White Balance : Manuel White Balance Electrical Shutter : 1/60 (PAL = 50) - 1/96000 Power Source With 21 mm Indian Coupler : DC 12 V 500 (SMPS)
Endo Camera (Digital Video Camera), ENT Equipment exporters, ENT Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)