

Product Name :
Defibrillator**Product Code :**
NMCBS139010**Description :**

Defibrillator Our inexplicable quality portable defibrillator is specially designed to automatically examine the ECG of heart patients. We are counted as the most preeminent external defibrillator manufacturer and one of the leading defibrillator exporters from India. A perfect combination of voice and text is in built to guide the user throughout the operation. It provides accurate results as LCD screen is incorporated in this. Light weight, versatile power supply and high portability are some of its unique features, which make it highly popular. Paramedic DEF-501 is intended for use during Advanced Cardiac Life Support training, Basic Life Support training, or any other physician approved trainings that deal with the operation of an automated external defibrillator. Paramedic DEF-502 is designed to automatically analyze the ECG of the patient. Analysis starts as soon as the ECG pads are properly attached. The user is guided by a combination of text and voice prompts throughout the rescue operation. The ECG of the patient is displayed in a liquid crystal display (LCD). Paramedic DEF-503, in automatic mode, is designed to automatically analyze the ECG of the patient. Analysis starts as soon as the ECG pads are properly attached. The user is guided by a combination of text and voice prompts throughout the rescue operation. The ECG of the patient is displayed in a liquid crystal display (LCD). Paramedic DEF-504 is a portable semi-automated external defibrillator with the following modes of operation : Automated External Defibrillation (AED) Mode Manual Mode AED and Manual Mode Defibrillation, Synchronized Cardio version, ECG Monitoring Mode (Lead III), SpO2 Monitoring, Heart Rate and SpO2 Alarm System, Lightweight and highly portable, Versatile power supply, Automatic and operator initiated self tests Paramedic DEF-505 defibrillator / monitor are designed to accommodate both basic and advanced life support personnel. AED and Manual Mode Defibrillation, Synchronized Cardio version, Reusable External Paddle (Adult, Pediatric) & Multifunction Defib. Pads ECG Monitoring Mode (3 Lead, 5 Lead), SpO2 Monitoring, Heart Rate and SpO2 Alarm System, Automatic and operator initiated self tests Defibrillator, ENT Equipment exporters, ENT Equipment manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)