

**Product Name :**  
Anti Bedsore Mattress

**Product Code :**  
HBMAT2010003



**Description :**

Anti Bedsore Mattress:  
Anti-Bedsore mattress (Air Jet Type).  
Tiny Holes on the surface for decreasing humidity.  
Replaces the foul air and reduce growth of microbes.  
It can effectively prevent from bedsore by dispersing weight of patient's body.  
Suitable for long term bedridden patients.

**Technical Specification :**

Very easy to operate - Adjustable comfort range with pressure-adjustable control knob  
Provides perfect support to bedridden patients  
Waterproof, antibacterial and fire resistant mattress offers maximum benefits  
It relieves the pain of sores and pressure ulcers by maintaining alternating pressure  
Perfect choice while you are looking for a portable, light weight and a perfect space saver option