

Product Name :
Auto Refractometer Matronix Q 30 Koria**Product Code :**
NMCBS139016**Description :**

Auto Refractometer Matronix Q 30 Koria We are engaged in offering Auto Refractometer Matronix q 30+ korean based technology. Technically advance with 3 year of warranty and money back garanty and has special IOL feature which enable to take the refractions of post operated patients. Features: $\dot{i}_{\dot{c}}\frac{1}{2}$ The most advanced technique of a graph edition, make the graph measured quickly and clear. $\dot{i}_{\dot{c}}\frac{1}{2}$ The precedent optical technique, consistent measurement and the accurate figure. Specifications: Refraction measurement Sphere (sph): -25D~+22D, 0.12D/0.25D steps Cylinder (cyl): -80D~ +8D, 0.12D/0.25D steps Axis : 0~180 $\dot{i}_{\dot{c}}\frac{1}{2}$, 1 $\dot{i}_{\dot{c}}\frac{1}{2}$ steps Pupil distance 45~85 mm, 0.5mm steps Vertical distance 0m, 12mm, 13.75mm Minimum pupil diameter 02.0mm Charter Automatic fogging method Display 5.7 inch true color LCD Print Import thermal printer Power Supply 100~120V/200~240V, 50/60HZ Dimension 275(W) * 509(D)*430-460(H)mm Weight 18Kg Auto Refractometer Matronix Q 30 Koria, ENT Equipment exporters, ENT Equipment manufacturers