

**Product Name :**  
A BRAND NEW DIRECT OPHTHALMOSCOPES

**Product Code :**  
NMCBS139027



**Description :**

A BRAND NEW DIRECT OPHTHALMOSCOPES OPHTHALMOSCOPES High quality optics for examining the interior structures of the eye, especially the retina, in a mini, pocket-size format OPHTHALMOSCOPES HEINE MINI 2000 LED critical tool offers all of the important features of a high-end ophthalmoscope combined in a compact, durable, portable configuration. Special Features: Enhanced XHL Xenon Halogen Technology; concentrated white light for perfect illumination & crisp retinal detail -20 to 20 with 18 lenses; -1 to +1 1/4 in single steps 5 apertures for full diagnostic capability H 6.75 in. x W 1.25 in. x D 0.9 in. Weight: 3.7 oz A BRAND NEW DIRECT OPHTHALMOSCOPES, ENT Equipment exporters, ENT Equipment manufacturers