

**Product Name :**  
TRIAL LENS SET**Product Code :**  
NMCBS139035**Description :**

TRIAL LENS SET We provide a range of Trial Lens Set that find applications in hospitals, private nursing home, clinics and other medical institutions. These lenses are available in different capacities that help the ophthalmologists in proper diagnosis and treatment of eye disease in patients. The trial lens case are made out of refined hard plastics that offers maximum protection to the lenses against impacts. Specifications: 32 Pairs Each cx and cc sph 20.00 Pairs Each cx and cc cyl 0.25 to 4.00 in 0.25 Steps 0.25 to 3.50 in 0.25 Steps 4.50 to 6.00 in 0.50 Steps 4.00 to 6.00 in 0.50 Steps 7.00 to 14.00 in 1.00 Steps Prisms 1/2, 1, 2, 3, 4, 5, 6, 8, 10, 12 16.00 to 20.00 in 2.00 Steps TRIAL LENS SET, ENT Equipment exporters, ENT Equipment manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,