

**Product Name :**  
ROTTER TEST TYPE DRUM

**Product Code :**  
NMCBS139037



**Description :**

ROTTER TEST TYPE DRUM Our range of Rotter Test Type Drum is has high strength and is available in wall type metal body. These are internally printed in accordance with snellen chart and our clients can avail in all the Indian languages. Various tests in which these can be availed are Duochrome test, Friend test, Spot light test and Worth four dots test. Offered at market leading prices, these meet the requirements of eye care hospitals and clinics. Features: Compact & light weight Remote control with Cord Pleasing color to match all interiors Include color deficiency test ROTTER TEST TYPE DRUM, ENT Equipment exporters, ENT Equipment manufacturers