

**Product Name :**  
DIGITAL PORTABLE VIEW TESTER

**Product Code :**  
NMCBS139039



**Description :**

DIGITAL PORTABLE VIEW TESTER The Digital Portable View Tester is distinguished by its high accuracy, wide measuring range, simplicity of operation & above all it is a really portable instrument operating on three pencil cells of size. we provide high quality tester to our clients. The measuring data of the digital refraction can also be printed by the instrument if the digital refraction has the data-output function. Features: 7.0 inch color LCD touch screen Stylish and compact design Auto detachment printer Independent PD adjustment Built-in LED for near vision illumination DIGITAL PORTABLE VIEW TESTER, ENT Equipment exporters, ENT Equipment manufacturers