



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Labour & Delivery Bed (With Accessories)

**Product Code :**  
NMCBS120019



**Description :**

Labour & Delivery Bed (With Accessories) Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes. Size : Body section: 1000L x 800W x 750H mm., Leg section: 900L x 800W x 750H mm. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width. Epoxy coated tubular steel frame. Carrying capacity: 150 - 180 kg. Supplied with S.S. bowl. Knockdown construction



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)