

Product Name :
Examination Couch Electric**Product Code :**
NMCBS119006**Description :**

Examination Couch Electric Electrically controlled, multi function, variable height couch has been designed to meet the ever growing needs of the busy practice. This new design has evolved from extensive user group feedback. Features Breathing hole option Swing Away siderails fitted as standard Variable height Electric hand-set control Adjustable head/back rest section 4 brake/swivel castors, for ease of manoeuvrability and safety Comfortable Upholster waterproof foam surface Accessories: couch end (dual use - push steer bar and patient safety feature), IV pole and complete lithotomy set (brackets, straps & pole) Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Piece Port of Dispatch: New Delhi Production Capacity: 500 units per month Delivery Time: As per buyer Packaging Details: Standard export packing