

**Product Name :**  
Examination Couch 3 Section**Product Code :**  
NMCBS119007**Description :**

Examination Couch 3 Section Height adjustment by electromotive drive  $\pm 1/2$  The electric motor is controlled by 1 foot switch (possible option: Hand switch)  $\pm 1/2$  Head part step less up- (+50 $\pm 1/2$ ) and downwards (-45 $\pm 1/2$ ) adjustable by gas spring.  $\pm 1/2$  Length head part: 450 mm  $\pm 1/2$  Middle part fixed. Length: 1.290 mm  $\pm 1/2$  Foot part step less adjustable from downwards vertical position (-90 $\pm 1/2$ ) to upwards position (+28 $\pm 1/2$ ) by gas spring.  $\pm 1/2$  Length foot part: 450 mm  $\pm 1/2$  Upholstery base: Laminated wood.  $\pm 1/2$  Length min.: 1.790 mm, Length max. 2.190 mm Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Piece Port of Dispatch: New Delhi Production Capacity: 100 units per month Delivery Time: As per buyer Packaging Details: Standard export packing