

**Product Name :**  
Blood Transfusion Chair**Product Code :**  
NMCBS122003**Description :**

Blood Transfusion Chair Blood sampling chair with fixed height, Height= 670 mm Comfortable upholstery with good density, width 550 Steel quality, epoxy paint Adjusting back rest via 1 gas jack Adjusting leg rest via a self locking mechanism system Equipped with 2 blood test splints, every position adjustable by serrated blocks system Load capacity: Dynamic 130kg Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 2 Pieces Port of Dispatch: New Delhi Production Capacity: 500 units per month Delivery Time: As per buyer Packaging Details: Standard export packing

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,