

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Blood Donor Chair

Product Code :
NMCBS122006



Description :

Blood Donor Chair Blood sampling chair with fixed height, Height= 460 mm Comfortable upholstery with good density, width 420 Steel quality, epoxy paint Adjusting back rest via 1 gas jack Equipped with 1 blood test splint, every position adjustable by serrated blocks system Load capacity: Dynamic 130kg Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 2 Pieces Port of Dispatch: New Delhi Production Capacity: 500 units per month Delivery Time: As per buyer Packaging Details: Standard export packing

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)