

**Product Name :**  
Blood Collection Chair**Product Code :**  
NMCBS122007**Description :**

Blood Collection Chair The frame is made of flat oval steel tubing The padded seat, back and armrests provide the patient with comfortable seating. The seat height is 500 mm. The distance between the wedge-shaped padded armrests 500 mm Maximum weight load 150 kg Knockdown construction Other Information Minimum Order Quantity: 2 Pieces Port of Dispatch: New Delhi Production Capacity: 500 units per buyer Delivery Time: As per buyer Packaging Details: Standard export packing Other Information Pay Mode Terms: L/C (Letter of Credit), D/P (Documents Against Payment), T/T (Bank Transfer) Minimum Order Quantity: 5 Pieces Port of Dispatch: New Delhi Production Capacity: 1000 units per month Delivery Time: As per buyer Packaging Details: Standard export packing