

**Product Name :**  
I.C.U. Multi-function Electric Bed With Weighing Scale**Product Code :**  
NMHOPB00007**Description :**

I.C.U. Multi-function Electric Bed With Weighing Scale Multi-Function ICU Bed electrically operated back rest tilting 0-65°, Foot rest tilting 0-43°, trendelenburg tilting 10°, Reverse trendelenburg tilting 10°, Height Adjustment (710mm to 530mm). Size: L 2330mm x W1050mm x H 530-710mm (Adjustable Height). All functions controlled with LINAK (Denmark) Power device (Actuator) In Build control panel on both side of Guard rail for Nurse/Patients to adjust bed position easily. Patient Lock-Out system for Nurse Staff at foot board. CPR position button for bed to reach Zero position quickly. Bed position data can be stored and recalled. Function of weighing scale, weight data can be stored and managed. Off-bed alarming. Removable High Quality ABS engineering Plastic head & foot panel. Size 1955mmX 900mm Four easy lifting Guard rail controlled by gas spring. Epoxy coated Mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 150mm dia swivel noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed