

**Product Name :**  
ICU Bed Electric Five Function**Product Code :**  
NMHOPB00009**Description :****ICU Bed Electric Five Function India China****Technical Specification :**

I.C.U. Bed, Electric, Five Function ICU Bed (Five function) with electrically operated back rest tilting 0-80

knee rest tilting 0-35

trendelenburg tilting 0- 20

Reverse trendelenburg tilting 0-20

Size: L 2140mm

W 940mm

H 590-990mm (Adjustable Height). All functions controlled with DEWERT/T-motion motor from Germany/Taiwan Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. Equipped with battery backup. 125mm dia noiseless castors with Central braking system. Back rest & Knee rest can be adjusted simultaneously with single button control. Provision for I.V. Rod on both sides of the bed

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE,](#)