

**Product Name :**  
Naugramedical ICU Bed Mechanical Five Function**Product Code :**  
NMHOPB00011**Description :**

Naugramedical ICU Bed Mechanical Five Function

**Technical Specification :**

Naugramedical ICU Bed Mechanical Five Function with mechanically operated back rest, knee rest tilting 0-35, trendelenburg 0-20/ reverse trendelenburg 0-20 (TR/RTR) and height position. Size: L 2140mm

W 980mm

H 420-720mm (Adjustable Height). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. C Type design. Epoxy coated mild steel frame work and 4 section perforated top. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed.