

Product Name :
ICU Beds Electric Three Function**Product Code :**
NMHOPB00014**Description :****ICU Beds Electric Three Function India China****Technical Specification :**

I.C.U. Bed, Electric, Three Function ICU Bed (3 function) with electrically operated remote controlled back rest 0-70

knee rest 0-35

and height position. Size: L 2160mm

W 940mm

H 520-770mm (Adjustable Height). All functions controlled with LINAK (Denmark) Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]