

**Product Name :**  
I.C.U. Bed, Electric, Three Function**Product Code :**  
NMHOPB00015**Description :**

I.C.U. Bed, Electric, Three Function Three-Function Electric Bed with electrically operated remote controlled back rest tilting 0-75°½, Knee rest tilting 0-35°½, Height Adjustment (640mm ÷ ½ 430mm). Size: L 2155mm ÷ ½ W 1050mm ÷ ½ H 430-640mm (Adjustable Height). All functions controlled with LINAK (Denmark) power device (Actuator) Removable High Quality ABS engineering Plastic head & foot panel, powder coated steel bed board with ventilated holes 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed