

Product Name :
ICU Bed 3 function Electric**Product Code :**
NMHOPB00016**Description :**

ICU Bed 3 function Electric with electrically operated remote controlled back rest 0-70°, knee rest 0-35° and height position. All functions controlled with LINAK (Denmark) Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Size: L 2140mm x 1/2 W 940mm x 1/2 H 500-750mm (Adjustable Height). Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.