

Product Name :
ICU Bed Electric Five Function with Weighing Scale System

Product Code :
NMHOPB00023



Description :

ICU Bed Electric Five Function with Weighing Scale System India China

Technical Specification :

I.C.U. Bed, Electric, Five Function (WITH WEIGHING SCALE SYSTEM) ICU Bed (Five function) with electrically operated back rest tilting 0-80

nee rest tilting 0-35

trendelenburg tilting 0-20

Reverse trendelenburg tilting 0-20

Size: L 2160mm

W 940mm

H 540-770mm (Adjustable Height). All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with battery backup & with Weighing scale system. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with bumpers protection. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,