

**Product Name :**  
Hospital ICU Bed Electric Five Function**Product Code :**  
NMHOPB00024**Description :****Hospital ICU Bed Electric Five Function India China****Technical Specification :**

I.C.U. Bed, Electric, Five Function ICU Bed (Five function) with electrically operated back rest tilting 0-80

knee rest tilting 0-35

trendelenburg tilting 0-20

Reverse trendelenburg tilting 0-20

Size: L 2160mm

W 940mm

H 480-730mm (Adjustable Height). All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with battery backup. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with bumpers protection. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE,](#)