

Product Name :
Medical ICU Bed Electric Five Function**Product Code :**
NMHOPB00020**Description :****Medical ICU Bed Electric Five Function India China****Technical Specification :**

I.C.U. Bed, Electric, Five Function ICU Bed (5 function) with electrically operated remote controlled back rest tilting 0-80

knee rest tilting 0-35

trendelenburg tilting 0-20

Reverse trendelenburg tilting 0-2

Size: L 2160mm

W 940mm

H 480-730mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,