

**Product Name :**  
ICU Bed Electric Four Function With Sitting Position**Product Code :**  
NMHOPB00025**Description :****ICU Bed Electric Four Function With Sitting Position India China****Technical Specification :**

I.C.U. Bed, Electric, Four Function With Sitting Position ICU Bed (4 function) with electrically operated remote controlled sitting position 0-70

back rest tilting 0-70

knee rest tilting 0-35

Size: L 2140mm

W 940mm

H 570-780mm (Adjustable Height). All functions controlled with electric Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. With Phenolic Laminating Sheet for X-ray. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section mesh mattress platform top. Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,