

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Anti Bedsore Mattress

**Product Code :**  
HBMAT2010009



**Description :**

Anti Bedsore Mattress Anti Bedsore Mattress (Air Jet Type). Tiny Holes on the surface for decreasing humidity. Replaces the foul air & reduce growth of microbes. Power supply: 220V/50Hz, 1A/250V, Air output: 3.5-4 liter/min. Pressure range: 50 mmHg to 90 mmHg. Material: PVC Size:200cm(L) x 90cm(W) x 6.5cm (H) Low Noisy, Long life time benefit

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)