

Product Name :
Stretcher Army 4 Fold With Telescopic Lifting Handles

Product Code :
NMSTR10024



Description :

Stretcher Army 4 Fold With Telescopic Lifting Handles Light Weight, easy clean, aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 220x54x15cm. Folded Size (LxWxH): 49x24x19cm. Self-weight: 11 kg, Load Bearing: 159 kg. Supplied with carry bag