

Product Name :
Salvage Stretcher Yellow**Product Code :**
NMSTR10026**Description :**

Salvage Stretcher Yellow Designed for fractured and immobile serious patients. Light Weight, easy clean, anodized aluminium alloy framework with safety belts. Size (LxWxH): Max. unfolded Size (LxWxH): 214x42x7cm, Min. unfolded Size (LxWxH): 167x42x7cm. Detachable and folding. This Salvage stretcher can be used with head immobilizer. Length of stretcher adjustable according to patient's body. The foot-end is designed as a narrow-frame structure. Self-weight: 8 kg, Load Bearing: 159 kg

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE