

Product Name :
Baby Bassinet**Product Code :**
NMINF10002**Description :**

Frame work made of ERW tubular pipe of 2.5 cms. dia. mounted on four 5 cms dia. castors.
Removable CRIB of Preplex transparent Crib with mattress and R-trend positions.
Trendlenburg / Rev. Trendlenburg.
Laminated cupboard for storage, I.V Rod.
Pre-treated and epoxy powder coated.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#).