

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Intravenous Transfusion Set

Product Code :
NMCBS35026



Description :

Intravenous Transfusion Set Extra strong, spike with/without air vent bacteria barrier at air inlet. Flexible drip chamber facilitate rapid adjustment of fluid level. Disc filter in drip chamber. Kink resistant non-toxic medical grade P.V.C. tubing. Strong roller clamp with a provision to hold the tube for changing the fluid container without disturbing the flow rate. Extra thick self sealing latex bulb. Luer lock connectors also available. E.T.O. Sterile and Pyrogen Free

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,